Instructions For Use

Table of Translations

English, English (US).....	3
20	Arabic, العربية
Български, Bulgarian	36
中文（简体）, Chinese Simplified.....	54
Hrvatski, Croatian	71
Čeština, Czech	86
Dansk, Danish	103
Nederlands, Dutch	120
Eesti keel, Estonian	140
Suomi, Finnish	155
Français, French	173
Deutsch, German	191
Ελληνικά, Greek	211
Magyar, Hungarian.....	227
Italiano, Italian.....	246
日本語, Japanese.....	263
한국어, Korean.....	278
Latviešu, Latvian.....	295
Lietuvių, Lithuanian	313
Polski, Polish	333
Português, Portuguese.....	349
Română, Romanian	367
Русский, Russian	386
Slovenský, Slovak.....	404
Slovenščina, Slovenian	421
Español, Spanish	440
Svenska, Swedish.....	458
Türkçe, Turkish	475

Instructions for Use

X-Ray Tube Assembly

English, English (U.S.)

Document Information	7
Legend	8
Intended Purpose	11
Intended Use	11
Intended User Population	11
Intended Context	11
Indications and Contraindications	11
Essential Performance	11
Medical Benefit	12
Safety Information	13
Safety Notice	13
Radiation Protection	14
Electrical Safety	14
Thermal Safety	14
Cooling / Insulation Oil	14
Electromagnetic Compatibility	15
Limits for Operation	15
Transportation and Storage	15
Disposal	16
Conformity	17
Installation and Maintenance	18
Installation	18

X-ray Tube Assembly Conditioning	18
General Maintenance	18
Corrective Maintenance	18
Planned Maintenance	19
Cleaning	19
Disinfection	19
Checks performed by System Operator	19
General Information	20
Carrier Damage	20
Warranty	20
Return Process	20
Service and Manufacturing Locations	21

1. Document Information

Document ID: 453577520510
Revision: B
Release Date: 2022-11

Copyright

© 2022 Chronos Imaging, LLC

All rights reserved. Reproduction in whole or in part is prohibited without the prior written consent of the copyright holder. Chronos Imaging, LLC, reserves the right to make changes in specifications or to discontinue any product at any time without notice or obligation, and is not liable for any consequences resulting from the use of this publication.

Equipment is subject to change without notice. All changes will be in compliance with regulations governing manufacture of medical equipment. Printed in the U.S.A. Document originally drafted in English.

2. Legend

Refer to Instruction Manual / Booklet

Warning: Ionizing Radiation

Warning: Electricity

Warning: Explosive Material

Warning: Non-Ionizing Radiation

General Warning Sign

Fragile, Handle with Care

Humidity Limitation

Atmospheric Pressure Limitation

Temperature Limit

Keep Dry


This Way Up


Stacking Limit by Mass

Protective Earth (Ground)

Earth (Ground)

Waste Electrical and Electronic Equipment (WEEE)

CE Mark

Intertek
NRTL Agency Approval

yyyy-mm

Address of Legal Manufacturer

yyyy: Year of Manufacture
mm: Month of Manufacture


Dimensions of Large Focal Spot

Dimensions of Small Focal Spot

Name of Product

Value of Permanent Filtration

Medical Device Symbol

Serial Number

Catalogue Number

Authorized Representative in the European Community

Unique Device Identifier

eIFU indicator

Consult Instructions for Use

X-Ray tube housing assembly

X-RAY TUBE

X-ray tube

Authorized Representative
in Switzerland

Nominal X-ray Tube Voltage

Value of nominal X-ray tube
voltage and respective
standard

3. Intended Purpose

Intended Use

This diagnostic X-ray tube housing assembly is an X-ray generating tube encased in a radiation-shielded housing that is intended for patient diagnostic imaging. The device is intended to be a replacement X-ray tube used in stationary installed C.T. Scanners that use the O.E.M. tubes for patient diagnostic imaging. The product does not come in direct contact with the patient or operator while in use. Only qualified field service technicians are authorized for installation, maintenance, and repair activities.

Intended User Population

Only qualified service technicians are authorized for installation, maintenance and repair activities.

Intended Context

This X-ray tube assembly is intended to be used in a non-condensing, climate controlled, indoor, clinical environment and is not intended to be used for homecare or by lay persons. The X-ray tube assembly is intended to be used in permanently installed stationary systems. The use frequency is not limited within the nominal expected lifetime. This X-ray tube assembly is intended to be re-used. The re-use of this X-ray tube assembly does not require special treatment. Additional treatment for re-use may be required and verified by the manufacturer of the X-ray system for specific applications.

Indications and Contraindications

Indications and contraindications are not defined at the X-ray tube housing assembly level. Indications and contraindications are defined by the manufacturer of the X-ray system in accordance with the intended use, the intended purpose and the medical purpose as described in the technical documentation of the X-ray system

Essential Performance

The X-ray tube assembly itself does not have essential performance, nor is any function of the X-ray tube assembly intended to contribute to the essential performance of the X-ray system it is installed on.

Medical Benefit

This X-ray tube assembly is a component part of a diagnostic X-ray imaging system. As an integral part of patient care, diagnostic X-ray is used in the diagnosis of serious diseases. Diagnostic X-ray is a safe, painless and cost-effective way of gathering information that may otherwise be unavailable or require a more expensive and/or riskier diagnostic test. It is particularly useful in emergency diagnosis and treatment. X-ray imaging is useful to diagnose bone injury and disease such as fractures, bone infections, arthritis and cancer. Diagnostic X-ray imaging procedure often identifies abnormalities early in the progression of a disease before some medical problems are apparent with other diagnostic tests. This early detection allows a disease to be treated earlier in its course. Diagnostic radiology involves some exposure to radiation. Special care is taken during the exam to ensure minimum exposure and maximum safety for the patient.

A medical procedure involving radiation should be done only when there is a justification. There should be an appropriate medical reason for the X-ray to be performed. The issue of medical radiation exposure is not a matter of safety. It is a matter of benefit compared with risk. That decision can only be made by someone who is familiar with the medical condition and the care that is necessary to properly manage it. The decision to have a medical test that involves radiation performed must be made collectively between the patient and his/her physician. All benefit/risk evaluations must be made on a case by case basis and are an inherent part of the medical art. Any diagnostic test should be justified by the risk of not having the test performed. This should be the basis for decisions made by physicians.

4. Safety Information

Safety Notice

As with any medical device, there are risks associated with use. Typical risks related to X-ray tube assemblies are X-ray tube arcing, oil leakage, use errors, biological hazards and energy hazards such as high voltage, low voltage, thermal and radiation hazards. These risks are inherent to X-ray tube assemblies and have been mitigated by conformity to harmonized standards. In addition, limiting use of the device to the useful lifetime, correct installation and verification of correct installation by the user ensure that any residual risks of the device are reduced as low as possible.

These Instructions for Use are designed to make it possible to work with the X-ray tube housing assembly in a safe manner. Operate the X-ray tube housing assembly only in compliance with the safety instructions in this manual and do not use it for purposes other than for which it is intended. The X-ray system may only be operated by qualified persons who have the necessary expertise in radiation protection and who have been instructed in how to operate the X-ray system.

It is always the user who is responsible for compliance with the regulations as applied to installation and operation of the X-ray unit

- You must never use the X-ray tube housing assembly or any associated diagnostic imaging equipment if it has any electrical, mechanical or radiological defects. This particularly applies to malfunctioning indicators, displays, warnings and alarms.
- Do not install the X-ray tube assembly to another equipment or medical device other than those with which it is compatible.
- Chronos Imaging, LLC is responsible for the safety features of its products only if it is maintained, repaired or modified by persons explicitly authorized to do so.
- As with any technical appliance, this equipment requires
 - correct operation
 - regular, competent maintenance
 - care
- If you operate the X-ray equipment, and hence the X-ray tube housing assembly, incorrectly or if the user fails to maintain it properly, Chronos Imaging, LLC cannot be held responsible for any malfunctions, damages or injuries.
- Safe operation of the X-ray tube assembly is only guaranteed when it is used according to its specification. In case the specification limits are disregarded, there is the danger of oil leaks and expelled component parts as a result of implosion/ explosion. In such cases the manufacturer of the X-ray tube assembly excludes any liability. Any guarantee claims for this product are rejected.
- The safety circuit of the X-ray tube assembly that prevents the switching ON of radiation when the specified temperature limit of the X-ray tube assembly is exceeded shall be connected before the X-ray tube assembly is set to work for the first time. It shall neither be removed nor modified.
- Any serious incident that occurs in relation to the X-ray assembly must be reported to Chronos Imaging, LLC and, if in the European Union, the competent authority of the Member State in which the user and/or patient is established.

Radiation Protection

Make certain before every X-ray exposure that all necessary radiation precautions have been taken, including confirming acceptable image quality.

You can find information about radiation precautions in the instructions for use for the individual X-ray systems with which you are using this X-ray tube assembly.

Electrical Safety

Only trained maintenance staff may remove the covers from the X-ray tube assembly.

This X-ray tube assembly may only be used in medical rooms which meet the requirements of relevant national and international standards and laws.

To avoid the risk of electric shock, this equipment must only be connected to a supply with protective earth.

Thermal Safety

To prevent damage due to thermal overloading, precautions must be taken to ensure that the X-ray tube assembly is not operated outside its specified load parameters. This prevention reduces risks to the patient, operating staff, third parties and the environment.

The x-ray tube housing assembly may be equipped with a thermal housing safety switch. Activation of the switch will inhibit additional exposures.

When this occurs, the unit will remain inoperable only until a safe operating temperature is restored, as determined by the imaging systems software. Depending on the cooling system, it may require 5 to 10 minutes of cooling time to reduce the temperature.

Cooling / Insulation Oil

If there is evidence of oil leakage from the X-ray tube housing assembly, heat exchanger, or hoses, immediately cease operation of the system and inform your service organization or supplier.

In the event of leakage or spillage of oil, ensure that the oil is removed as soon as the equipment has ceased operation and the area is safe. Use liquid-absorbent material and dispose in accordance with local environmental laws and regulations.

Electromagnetic Compatibility (E.M.C.)

In accordance with its intended use, this electronic component designed as part of a diagnostic imaging system according to the laws governing E.M.C., which defines the permitted emission levels from electronic equipment and its required immunity against electromagnetic fields.

Electronic apparatus that satisfies the E.M.C. requirements is designed so that under normal conditions there is no risk of malfunction caused by electromagnetic interference. However, in the case of radio signals from high-frequency transmitters with a relatively high transmitting power, the risk of electromagnetic incompatibility when operated in close proximity to electronic apparatus cannot be totally ruled out.

In unusual circumstances unintended functions of the apparatus could be initiated, possibly giving rise to undesirable risks for the patient or user. For this reason, all kinds of transmission with mobile radio equipment should be avoided. This also applies when the apparatus is in STANDBY mode. Mobile telephones must be switched OFF in designated problem zones.

The responsibility to guarantee compliance with the electromagnetic emission limit values and the fulfillment of all associated regulations and laws lies with the company responsible for the installation of the X-ray tube assembly.

Limits for Operation

Safe operation depends on ensuring the X-ray tube assembly is used within its limits. Failure to obey these limits can result in the danger of oil leaks. This loss of cooling media can cause the X-ray tube assembly to overheat which can lead to expelled component parts as a result of implosion / explosion. Refer to the Environmental Limits for Operation in the supplemental I.F.U. documents.

- Do not operate the X-ray tube assembly in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.
- Do not use detergents or disinfectants, including those used on the patient, that can create explosive gas mixtures.

Transportation and Storage

Safe operation depends on ensuring the X-ray tube assembly is used within its limits. Failure to obey these limits can result in the danger of oil leaks. This loss of cooling media can cause the X-ray tube assembly to overheat which can lead to expelled component parts as a result of implosion/ explosion.

It is recommended to transport x-ray tube housing assemblies in the original packaging material. Transportation using other packaging material or in any other way may result in injury and/or product damage and may void the warranty.

Refer to the Environmental Limits for Transportation detailed in the supplementary I.F.U. documents.

Disposal

The take-back obligation, correct disposal and recovery of the X-ray tube assembly refer to the European Waste Electrical and Electronic Equipment (WEEE) directive. They also refer to the requirements of local and transregional legal requirements.

The X-ray tube assembly manufacturer assembles state-of-the-art X-ray tube assemblies in terms of safety and environmental protection. If no parts of the X-ray tube assembly are opened and if the X-ray tube assembly is used correctly there are no risks to persons or the environment.

To obey regulations, sometimes it is necessary to use materials that are harmful to the environment. Discard these materials in a correct manner.

This X-ray tube assembly contains materials that are toxic. Do not discard the X-ray tube assembly together with industrial or domestic waste.

The manufacturer

- supports you in the disposal of the X-ray tube assembly in accordance with valid legal requirements.
- takes back the X-ray tube assembly.
- returns re-usable parts to the production cycle. Extensive test and quality assurance procedures as well as detailed checks of the components qualify these parts to meet the same high level of quality and functionality that is expected from new materials.
- makes a contribution to the protection of the environment.

In case you have questions concerning safe disposal, please consult the manufacturer in full confidence.

5. Conformity

The X-ray tube housing assembly meets the provisions of European CE Marking of Conformity bears the CE Mark Label. The assembly meets the provisions of F.D.A. 21 CFR 1020.30.

Refer to the Dimensional Data and Labeling Section for the location of C.E. Label on the X-ray Tube Housing Assembly.

CE 2797

UK
CA 0086

6. Installation and Maintenance

Installation

This X-ray tube is designed, manufactured, and validated to be compatible with specific diagnostic imaging equipment. Prior to installation, consult the 3rd Party Compatibility Table in the supplemental I.F.U. documents.

Only trained and qualified service personnel must perform installation of the X-ray tube housing assembly. The installation must be performed in accordance with the system documentation provided by the Original Equipment Manufacturer's installation, test and operating procedures. Consult the appropriate documentation for the removal and installation procedure.

This X-ray tube housing assembly is intended to be installed on a diagnostic imaging system that is fitted with an appropriate audible warnings and/or visual displays, where practicable, while radiation is being emitted.

Failure to follow the O.E.M. procedure for removal and installation of the X-ray tube housing assembly may result in injury to user or patient, damage to equipment and will void the warranty.

X-ray Tube Assembly Conditioning

It is imperative for proper operation and extended life that the X-ray tube housing assembly be properly conditioned after a period of non-use. This procedure is defined in the Original Equipment Manufacturer's diagnostic imaging system operator's manual. Failure to follow this procedure may reduce the life of the X-ray tube housing assembly and may void the warranty terms.

General Maintenance

The responsibility to recommend preventive maintenance in accordance with all applicable regulations and laws lies with the company responsible for the installation of the X-ray tube assembly.

As with any technical appliance these X-ray tube assemblies also require:

- Regular checks by the operator
- Regular planned and corrective maintenance

Corrective Maintenance

Corrective maintenance on X-ray tube assemblies is allowed to be performed only by the manufacturer of the X-ray tube assembly.

Faulty components of the X-ray equipment which affect the safety of the X-ray tube assembly must be replaced by genuine spare parts.

Planned Maintenance – General

X-ray tube assemblies contain mechanical components, which are subject to normal wear due to operation.

The correct setting of the electromechanical and electronic assemblies safeguards the functioning, image quality, electrical safety and radiation exposure of the patients, operating staff, third parties and the environment.

By taking these precautions, you maintain the operability and operational reliability of the system. As the user of the X-ray unit, you are obliged according to accident prevention regulations, the medical products law and other regulations to perform such precautionary actions.

Maintenance consists of tests that the user of the diagnostic imaging system can perform and maintenance that is performed under service agreements or by persons explicitly authorized to do so.

The user must check the X-ray equipment for apparent defects. If operational defects or other departures from normal operational behavior occur, the X-ray unit must be switched off and the appropriate service organization informed. Operation of the X-ray equipment may only be resumed when repairs are completed. Operation using faulty components may lead to an increased safety risk or unnecessarily high exposure to radiation.

It is recommended that you perform the tests indicated in the table on a regular basis and have the unit serviced by authorized representatives at least once a year. For heavily used equipment, preventative maintenance should be scheduled more often.

These precautionary measures prevent personal injury to patient and operator.

Cleaning

Cleaning of the X-ray tube assembly is not necessary as it is behind the covers of the X-ray imaging system. Cleaning of the X-ray tube assembly should not be performed.

Disinfection

Disinfection of the X-ray tube assembly is not necessary as it is behind the covers of the X-ray imaging system, unless otherwise required by the manufacturer of the X-ray imaging system.

Checks performed by System Operator

Interval	Scope of Work
Each use	Check of system error messages Check for visible oil leaks or other contamination
Daily	Check for damaged parts, missing labels and warning plates
Weekly	Check of all cables and connections (loose, damaged or broken)
Weekly	Check for unusual noises
Per system manual and/or relevant local or transregional standards and laws	Stability test
Per system manual and/or relevant local or transregional standards and laws	Check of image quality

7. General Information

Carrier Damage

In the unlikely event of damage due to handling by the carrier, it is important to follow these steps in order to receive proper credit:

1. Inspect the X-ray tube housing assembly immediately upon receipt. Check both the packing and the product for physical damage.
2. If there is physical damage, call the carrier immediately and order a “Joint Inspection” of both the packing and the product.

Within the United States:

- If the X-ray tube housing assembly has been shipped freight pre-paid, return the product in accordance with the instruction on the Return Form along with a copy of the Joint Inspection Report to Chronos Imaging, LLC.
- If the X-ray tube housing assembly has been shipped via the customer's carrier, the customer must initiate the claim process directly with the carrier.

Outside the United States:

- Please consult your local sales and service office for instruction for return.

Warranty

Warranties vary according to the specific X-ray tube assembly model. Contact your local representative for detailed information. Written copies are available upon request.

Return Process

Use the re-usable packaging of the new delivered X-ray tube assembly for the return shipment of the decommissioned X-ray tube assembly. This guarantees a transport of the returning X-ray tube assembly in packaging designed for this purpose.

Follow the instructions on the Return Form, fill in all required information and send it together with the X-ray tube assembly being returned.

Outside the United States:

If returning a product from outside the United States, please contact your local sales or service representative for return instructions.

8. Service and Manufacturing Locations

Device Manufacturer

Chronos Imaging, LLC
555 North Commerce Street
Aurora, IL 60586 USA

Website: www.chronosimaging.com

Authorized Representative in the European Community

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Authorized Representative in Switzerland

MDSS CH GmbH
Laurenzenvorstadt 61
5000 Aarau, Switzerland

UK Responsible Person

MDSS-UK RP Ltd.
6 Wilmslow Road
Manchester M14 5TP
United Kingdom

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All rights reserved. 저작권 소유자의 사전 서면 동의 없이 전체 또는 부분을 복제하는 행위는 금지됩니다. Chronos Imaging, LLC 는 사전 통지 또는 의무 없이 언제든지 사양을 변경하거나 제품 생산을 중단할 수 있으며, 본 문서를 사용함으로써 발생하는 결과에 대해 책임을 지지 않습니다.
MDSS GmbH